



Alcohol Screening Tests



How it works

The Breath Alcohol Detector is for rapid detection of the presence of alcohol in exhaled breath and provides relative Blood Alcohol Concentration (BAC) at 0.05%. The test is based on a chromogenic reaction. Alcohol, if present in the exhaled breath, reacts with the chemically coated crystals and produces a colour change. This colour change is proportional to the concentration of alcohol in the breath, which is an approximation of relative Blood Alcohol Concentration (BAC).

Contents

Instruction Sheet	1
Test per pack	2
Effectiveness	0.05%

Packing Details

Product Code	Description	Tests per Retail Pack	Retail Packs per Inner	Retail Packs per Shipper
SALCHO 104-2	SURESIGN Alcohol Screening	2	6	288



Sales - UK
The Miles Group Ltd.
Tel 01484 536344
Fax 01484 536220
info@milesgroup.co.uk
www.milesgroup.co.uk

Sales - Worldwide
Ciga Healthcare Ltd
Tel 02825 685385
Fax 02825 685145
info@cigahealthcare.com
www.cigahealthcare.com



Family Health Screening Tests



How it works

Family Health Screening Test strips are plastic strips onto which six separate reagent pads are affixed. The strips are intended for testing urine in the home for Nitrate, Protein, Urobilinogen, Blood, Specific Gravity and Bilirubin. The tests can be used to give a general indication of health and helps in the detection of diseases and infections that can effect the kidneys, urinary tract, liver & heart.

Contents

Instruction Sheet	1
Colour Chart	1
Test per Pack	6

Packing Details

Product Code	Description	Tests per Retail Pack	Retail Packs per Inner	Retail Packs per Shipper
SFHSC 101	SURESIGN Family Health Screening	6	6	288



Sales - UK
 The Miles Group Ltd.
 Tel 01484 536344
 Fax 01484 536220
 info@milesgroup.co.uk
 www.milesgroup.co.uk

Sales - Worldwide
 Ciga Healthcare Ltd
 Tel 02825 685385
 Fax 02825 685145
 info@cigahealthcare.com
 www.cigahealthcare.com